



Moticon SCIENCE

Price List

Effective July 19, 2019

Content

- 03 Sensor Insoles
Spare Parts
- 04 Power Supply
Product Bags
- 05 Desktop Software
Mobile App
- 06 Services
Packaging & Freight
- 07 Imprint

Sensor Insoles

Product	Description	PU	Price (Euro)
SCIENCE Sensor Insole	<p>Connectivity Bluetooth LE 4.2 or higher</p> <p>Sample Rates 10, 25, 50, 100 Hz</p> <p>Sensors 16 pressure sensors 6-axis IMU</p> <p>Memory 16MB internal</p> <p>Data Acquisition Recording, SmartRecording, Live, SmartSleep</p> <p>Data Output Pressure, total force, COP, acceleration, angular rates, time</p> <p>Power Supply Li-ion 2032 coin cell battery, 3.7V</p> <p>Sizes 9 double sizes (EU 32/33 - 48/49)</p>	1 pr.	1495.00
SCIENCE Dummy Insole	<p>Without integrated measurement technology</p> <p>Can be used as placebo for control group in intervention studies</p> <p>Can be used for testing shape and fit</p>	1 pr.	195.00

Spare Parts

Product	Description	PU	Price (Euro)
Battery Lid	Replacement, fits SCIENCE Sensor Insoles	1 pr.	20.00

Power Supply

Product	Description	PU	Price (Euro)
Coin Cell Battery	PD2032, Li-ion, rechargeable Voltage 3.7V, capacity 75mAh	1 pr.	20.00
Coin Cell Charger	2 charging slots for PD2032 coin cells Charge voltage 4.2V, Supply voltage 5.0V Micro-USB Dimensions 46x34x23 mm LED charge status indicator	1 pc.	95.00

Product Bags

Product	Description	PU	Price (Euro)
Product Bag	Compartments for 6 pairs of sensor insoles 2 flat compartments, 1 large side compartment Outer material Cordura, carry strap, zippers Dimensions 430x300x60 mm	1 pc.	145.00
Coin Cell Bag	Compartments for 20 PD2032 coin cells Outer material Cordura, zippers Dimensions 150 x 90 x 15 mm	1 pc.	25.00
Sensor Insole Bag	Fits 1 pair of sensor insoles Outer material mesh textile Dimensions 310 x 120 mm	1 pc.	15.00

Desktop Software

Product	Description	PU	Price (Euro)
SCIENCE Desktop Software ¹	Basic data management & analysis	1 lcs.	995.00
SCIENCE Analyze Module ²	Advanced data labeling, analysis & processing	1 lcs.	995.00
SCIENCE Video Module ²	Automated video synchronization Multi video channel management	1 lcs.	995.00
SCIENCE Gait Report ²	Automated gait report generation Custom report setup & export functions	1 lcs.	495.00
SCIENCE Jump Report ²	Automated jump report generation Custom report setup & export functions	1 lcs.	495.00
SCIENCE Balance Report ²	Automated balance report generation Custom report setup & export functions	1 lcs.	495.00

¹ Minimum system requirement: Windows 10, Intel i5 or equivalent, 4 GB RAM, 1 GB HDD, DirectX9

² Requires SCIENCE Desktop Software license

Mobile App

Product	Description	PU	Price (Euro)
SCIENCE Mobile App ³	Sensor insole control and settings Measurement organization, data transfer to Moticon Cloud or SCIENCE Desktop App	1 lcs.	0.00

³ For Android mobile phones, available in Google PlayStore

Minimum system requirements: Android 7 (Nougat), BLE 4.1, recommended system requirements: Android 9 (Pie), BLE 5.0

Services

Product	Description	PU	Price (Euro)
Support Service	12 months customer support Mandatory for 1st order	1 pc.	295.00
Software Update Service	12 months updates for SCIENCE Desktop App	1 pc.	295.00
Cloud Storage Service	12 months cloud storage, unlimited measurements	1 lcs.	295.00
Moticon Sensor Insole Replacement Service	On-demand service linked to purchase of each pair of Moticon Sensor Insoles Valid 12 months after purchase, unlimited exchanges	1 pr.	1195.00

Packaging & Freight

Product	Description	PU	Price (Euro)
Packaging ¹	Packaging fee	1 pc.	15.00
Freight National ²	Carriage fee Germany Carrier DHL Express	1 pc.	15.00
Freight Europe ²	Carriage fee for countries of European Union Carrier DHL Express	1 pc.	45.00
Freight International ²	Carriage fee all countries outside European Union Carrier DHL Express Includes custom declaration	1 pc.	150.00

¹ Per package, large orders may require multiple packages

² Per shipment, large orders may require multiple shipments

Imprint

Legal Notice

This document is the property of Moticon ReGo AG.
Any reproduction of this document in part or in whole is strictly prohibited.
For educational purposes only.
Errors and omissions excepted and subject to change without notice.

Notes on Prices

This price list informs you about delivery ex works (EXW).
Prices are net plus VAT.
For international orders, the VAT must be paid locally by the customer.
For non-EU orders, customs duty fees may apply according to national import tax rates.
The Moticon general terms and conditions (GTC) under www.moticon.de apply to this price list.
Moticon offers international carriage, insurance and custom services, but cannot be held liable for any loss that may occur during the carriage.

Further Information

For more detailed product information, please visit
www.moticon.de/insole3-overview
www.moticon.de/apps
www.moticon.de/accessories

Release Notes

Version 03.00.03
Date July 19, 2019

Vendor Information

Moticon ReGo AG
Machtlfinger Str. 21
81379 Munich
Germany

Phone +49 89 2000 301 60
EMail sales@moticon.de
Web www.moticon.de



E-Mail sales@moticon.de
Web www.moticon.de

